**AGITADOR ROTATORIO
ROLL & TILT MIXER
AGITATEUR ROTATOIRE**

Modelo / Model / Modèle 683

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Thank you for choosing a roll and tilt Nahita mixer, its high mechanical quality will allow you to have this equipment available for many years. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

In order to use this instrument, it is advisable to follow the following points:

- ◆ This manual should be available for all the users of this equipment.
- ◆ Any doubt can be clarified getting into touch with the NAHITA Technical Service.
- ◆ As an accuracy instrument, it should be used carefully avoiding the sudden movements and knocks.
- ◆ Never dismantle the different pieces of the magnetic stirrer to repair them, since it could produce a defective use of the whole equipment and a loss of the product warranty.

ATTENTION!!

IF AN EQUIPMENT IS NOT PROPERLY CLEANED AND DISINFECTED IT WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish.....	1-3
English	4-5
French.....	6-7

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	5
2. TECHNICAL SPECIFICATIONS	5
3. OPERATING INSTRUCTIONS	5

1. USES OF THE INSTRUMENT

The roll and tilt Nahita mixers are powerful and silent equipments based on a circular movement and rocking of 7 rotary rollers. With this equipment we can achieve a comfortable and quick way of levelling the disolutions of samples.

2. TECHNICAL SPECIFICATIONS

<i>Speed of rocking and rotation</i>	72 r.p.m
<i>Length of the roller</i>	240 mm
<i>Diameter of the roller</i>	35 mm
<i>Power Consume</i>	12 W
<i>Feeding</i>	220 V 50 Hz
<i>Dimensions</i>	395 x 250 x 100 mm
<i>Weight</i>	4 Kg

3. OPERATING INSTRUCTIONS

- ◆ Please unwrap the Nahita roll and tilt mixer avoiding sudden movements and protecting it from bumps, please put the equipment on a plane stable surface.
- ◆ The roll and tilt mixer is supplied with a Schuko standard wire this should be plugged to an earth wire. This should be handy so that all users can unplug the equipment when necessary.

The manufacturer and the distributor will not assume any responsibility for the damages produced to the equipment in its installing or damages to persons suffered by the wrong use of the electric connection. The tension should be 220V ±10% 50-60Hz.

- ◆ Once the equipment is connected a red pilot light will go on this pilot is situated on the superior surface of the equipment, it indicates that the equipment is connect to the electric net, to turn it on press the On/Off switch that is on the superior surface, with the equipment functioning the pilot will keep with the light on.
- ◆ The equipment has a protective fuse, if you need to substitute it, unscrew the top side of the device that carries the fuses and substitute the blown fuse.

Note: Take special care in not spilling liquid on the electric components of the equipment because it can suffer damages. If accidentally you spill liquid dry it as soon as possible this way you will avoid future problems.